## Package: pcv (via r-universe)

August 30, 2024

Version 1.1.0

Date 2023-08-12

Title Procrustes Cross-Validation

Maintainer Sergey Kucheryavskiy <svkucheryavski@gmail.com>

**Description** Implements Procrustes cross-validation method for Principal Component Analysis, Principal Component Regression and Partial Least Squares regression models. S. Kucheryavskiy (2023) <doi:10.1016/j.aca.2023.341096>.

**Encoding** UTF-8

License MIT + file LICENSE

Imports graphics, grDevices, stats

RoxygenNote 7.2.3

Suggests testthat

NeedsCompilation no

**Depends** R (>= 3.5.0)

URL https://github.com/svkucheryavski/pcv

BugReports https://github.com/svkucheryavski/pcv/issues

Repository https://svkucheryavski.r-universe.dev

RemoteUrl https://github.com/svkucheryavski/pcv

RemoteRef HEAD

RemoteSha 41c4399a499312a85d184b2a6c58a8c772899435

## Contents

corn	•	•			•	 		•			•	•	 •	 		•	•	•	•		•		•	•	•	2
getCrossvalParams	•		•	•	•	 			•	•	•	•		 			•			•				•		3
getR	•		•	•	•	 			•	•	•	•		 			•			•				•		3
getxpvorth	•		•	•	•	 		•	•	•	•	•	 	 			•	•		•				•		4
normalize	•		•	•	•	 		•	•	•	•	•	 	 			•	•		•				•		4

#### corn

		15
simpls	•••	14
rotationMatrixToX1		
plotD		12
pcvreg		12
pcvpls		10
pcvpcr		8
рсурса		6
peverossval		5
pcv		5

## Index

corn

Corn data

## Description

NIR spectra and moisture of 80 Corn samples.

#### Usage

data(corn)

## Format

A list with two matrices, spectra and moisture.

## Details

This is a part of Corn dataset, which was downloaded from Eigenvector Research, Inc. website (https://eigenvector.com/resources/data-sets/), where it is available publicly. This dataset contains several NIR spectra of corn samples recorded using different instruments. For our examples we took "mp5" spectra and corrected them using Standard Normal Variate transformation.

#### Source

1. Eigenvector Research, Inc. (https://eigenvector.com/resources/data-sets/)

getCrossvalParams Returns parameters for cross-validation based on 'cv' value

## Description

Returns parameters for cross-validation based on 'cv' value

## Usage

getCrossvalParams(cv, nobj)

## Arguments

CV	settings for cross-validation provided by user
nobj	number of objects in calibration set

getR	Creates rotation matrix to map a set vectors base1 to a set of vectors
	base2.

## Description

In both sets vectors should be orthonormal.

## Usage

getR(base1, base2)

## Arguments

base1	Matrix (JxA) with A orthonormal vectors as columns to be rotated (A $\leq J$ )
base2	Matrix (JxA) with A orthonormal vectors as columns, base1 should be aligned with

## Value

Rotation matrix (JxJ)

getxpvorth

## Description

Generates the orthogonal part for Xpv

## Usage

getxpvorth(q.k, X.k, PRM)

## Arguments

q.k	vector with orthogonal distances for cross-validation set for segment k
X.k	matrix with local validation set for segment k
PRM	projecton matrix for orthogonalization of residuals

## Value

A matrix with orthogonal part for Xpv

normalize	Normalization rows or columns of a matrix	
-----------	---	--

## Description

Normalization rows or columns of a matrix

## Usage

```
normalize(
   X,
   dim = 1,
   weights = if (dim == 1) 1/sqrt(rowSums(X^2)) else 1/sqrt(colSums(X^2))
)
```

## Arguments

Х	matrix with numeric values
dim	which dimension to normalize (1 for rows, 2 for columns)
weights	vector with normalization weights, by default 2-norm is used

5

## Description

Compute matrix with pseudo-validation set

#### Usage

```
pcv(
  X,
  ncomp = min(round(nrow(X)/nseg) - 1, col(X), 20),
  nseg = 4,
  scale = FALSE
)
```

#### Arguments

Х	matrix with calibration set (IxJ)
ncomp	number of components for PCA decomposition
nseg	number of segments in cross-validation
scale	logical, standardize columns of X prior to decompositon or not

## Details

This is the old (original) version of PCV algorithm for PCA models. Use pcvpca instead. Ane check project web-site for details: https://github.com/svkucheryavski/pcv

The method computes pseudo-validation matrix Xpv, based on PCA decomposition of calibration set X and systematic (venetian blinds) cross-validation. It is assumed that data rows are ordered correctly, so systematic cross-validation can be applied

#### Value

```
Pseudo-validation matrix (IxJ)
```

pcvcrossval

Generate sequence of indices for cross-validation

## Description

Generates and returns sequence of object indices for each segment in random segmented cross-validation

рсv

#### Usage

pcvcrossval(cv = 1, nobj = NULL, resp = NULL)

#### Arguments

CV	cross-validation settings, can be a number, a list or a vector with integers.
nobj	number of objects in a dataset
resp	vector or matrix with response values to use in case of venetian blinds

#### Details

Parameter 'cv' defines how to split the rows of the training set. The split is similar to cross-validation splits, as PCV is based on cross-validation. This parameter can have the following values:

1. A number (e.g. 'cv = 4'). In this case this number specifies number of segments for random splits, except 'cv = 1' which is a special case for leave-one-out (full cross-validation).

2. A list with 2 values: 'list("name", nseg)'. In this case '"name"' defines the way to make the split, you can select one of the following: '"loo"' for leave-one-out, '"rand"' for random splits or '"ven"' for Venetian blinds (systematic) splits. The second parameter, 'nseg', is a number of segments for splitting the rows into. For example, 'cv = list("ven", 4)', shown in the code examples above, tells PCV to use Venetian blinds splits with 4 segments.

3. A vector with integer numbers, e.g. 'cv = c(1, 2, 3, 1, 2, 3, 1, 2, 3)'. In this case number of values in this vector must be the same as number of rows in the training set. The values specify which segment a particular row will belong to. In case of the example shown here, it is assumed that you have 9 rows in the calibration set, which will be split into 3 segments. The first segment will consist of measurements from rows 1, 4 and 7.

#### Value

vector with object indices for each segment

рсvрса

Procrustes cross-validation for PCA models

#### Description

Procrustes cross-validation for PCA models

#### Usage

```
pcvpca(
   X,
   ncomp = min(nrow(X) - 1, ncol(X), 30),
   cv = list("ven", 4),
   center = TRUE,
   scale = FALSE,
   cv.scope = "global"
)
```

#### pcvpca

#### Arguments

Х	matrix with predictors from the training set.
ncomp	number of components to use (more than the expected optimal number).
CV	which split method to use for cross-validation (see description for details).
center	logical, to center or not the data sets
scale	logical, to scale or not the data sets
cv.scope	scope for center/scale operations inside CV loop: 'global' — using globally computed mean and std or 'local' — recompute new for each local calibration set.

#### Details

The method computes Procrustes validation set (PV-set), matrix Xpv, based on PCA decomposition of calibration set 'X' and cross-validation. See description of the method in [1].

Parameter 'cv' defines how to split the rows of the training set. The split is similar to crossvalidation splits, as PCV is based on cross-validation. This parameter can have the following values:

1. A number (e.g. 'cv = 4'). In this case this number specifies number of segments for random splits, except 'cv = 1' which is a special case for leave-one-out (full cross-validation).

2. A list with 2 values: 'list("name", nseg)'. In this case '"name"' defines the way to make the split, you can select one of the following: '"loo"' for leave-one-out, '"rand"' for random splits or '"ven"' for Venetian blinds (systematic) splits. The second parameter, 'nseg', is a number of segments for splitting the rows into. For example, 'cv = list("ven", 4)', shown in the code examples above, tells PCV to use Venetian blinds splits with 4 segments.

3. A vector with integer numbers, e.g. 'cv = c(1, 2, 3, 1, 2, 3)'. In this case number of values in this vector must be the same as number of rows in the training set. The values specify which segment a particular row will belong to. In case of the example shown here, it is assumed that you have 9 rows in the calibration set, which will be split into 3 segments. The first segment will consist of measurements from rows 1, 4 and 7.

Parameter 'cv.scope' influences how the Procrustean rule is met. In case of "global" scope, the rule will be met strictly - distances for PV-set and the global model will be identical to the distances from conventional cross-validation. In case of "local" scope, every local model will have its own center and scaling factor and hence the rule will be almost met (the distances will be close but not identical).

#### Value

Matrix with PV-set (same size as X)

#### References

1. S. Kucheryavskiy, O. Rodionova, A. Pomerantsev. Procrustes cross-validation of multivariate regression models. Analytica Chimica Acta, 1255 (2022) [https://doi.org/10.1016/j.aca.2023.341096]

#### Examples

```
# load NIR spectra of Corn samples
data(corn)
X <- corn$spectra
# generate Xpv set based on PCA decomposition with A = 20 and venetian blinds split with 4 segments
Xpv <- pcvpca(X, ncomp = 20, center = TRUE, scale = FALSE, cv = list("ven", 4))
# show the original spectra and the PV-set (as is and mean centered)
oldpar <- par(mfrow = c(2, 2))
matplot(t(X), type = "1", lty = 1, main = "Original data")
matplot(t(Xpv), type = "1", lty = 1, main = "PV-set")
matplot(t(scale(X, scale = FALSE)), type = "1", lty = 1, main = "Original data (mean centered)")
matplot(t(scale(Xpv, scale = FALSE)), type = "1", lty = 1, main = "PV-set (mean centered)")
par(oldpar)
```

```
pcvpcr
```

Procrustes cross-validation for PCR models

#### Description

Procrustes cross-validation for PCR models

#### Usage

```
pcvpcr(
  X,
  Y,
  ncomp = min(nrow(X) - 1, ncol(X), 30),
  cv = list("ven", 4),
  center = TRUE,
  scale = FALSE,
  cv.scope = "global"
)
```

#### Arguments

Х	matrix with predictors from the training set.
Υ	vector with response values from the training set.
ncomp	number of components to use (more than the expected optimal number).
cv	which split method to use for cross-validation (see description of method 'pcv-pls()' for details).
center	logical, to center or not the data sets
scale	logical, to scale or not the data sets
cv.scope	scope for center/scale operations inside CV loop: 'global' — using globally computed mean and std or 'local' — recompute new for each local calibration set.

8

#### pcvpcr

#### Details

The method computes pseudo-validation matrix Xpv, based on PCR decomposition of calibration set X, y and cross-validation.

Parameter 'cv' defines how to split the rows of the training set. The split is similar to cross-validation splits, as PCV is based on cross-validation. This parameter can have the following values:

1. A number (e.g. 'cv = 4'). In this case this number specifies number of segments for random splits, except 'cv = 1' which is a special case for leave-one-out (full cross-validation).

2. A list with 2 values: 'list("name", nseg)'. In this case '"name"' defines the way to make the split, you can select one of the following: '"loo"' for leave-one-out, '"rand"' for random splits or '"ven"' for Venetian blinds (systematic) splits. The second parameter, 'nseg', is a number of segments for splitting the rows into. For example, 'cv = list("ven", 4)', shown in the code examples above, tells PCV to use Venetian blinds splits with 4 segments.

3. A vector with integer numbers, e.g. 'cv = c(1, 2, 3, 1, 2, 3, 1, 2, 3)'. In this case number of values in this vector must be the same as number of rows in the training set. The values specify which segment a particular row will belong to. In case of the example shown here, it is assumed that you have 9 rows in the calibration set, which will be split into 3 segments. The first segment will consist of measurements from rows 1, 4 and 7.

Parameter 'cv.scope' influences how the Procrustean rule is met. In case of "global" scope, the rule will be met strictly - error of predictions for PV-set and the global model will be identical to the error from conventional cross-validation. In case of "local" scope, every local model will have its own center and hence the rule will be almost met (the errors will be close but not identical).

## Value

Pseudo-validation matrix (same size as X) with an additional attribute, 'D' which contains the scaling coefficients (ck/c)

#### References

1. S. Kucheryavskiy, O. Rodionova, A. Pomerantsev. Procrustes cross-validation of multivariate regression models. Analytica Chimica Acta, 1255 (2022) [https://doi.org/10.1016/j.aca.2023.341096]

#### Examples

```
# load NIR spectra of Corn samples
data(corn)
X <- corn$spectra
Y <- corn$moisture</pre>
```

```
# generate Xpv set based on PCA decomposition with A = 20 and venetian blinds split with 4 segments
Xpv <- pcvpcr(X, Y, ncomp = 20, center = TRUE, scale = FALSE, cv = list("ven", 4))</pre>
```

```
# show the original spectra and the PV-set (as is and mean centered)
oldpar <- par(mfrow = c(2, 2))
matplot(t(X), type = "1", lty = 1, main = "Original data")
matplot(t(Xpv), type = "1", lty = 1, main = "PV-set")
matplot(t(scale(X, scale = FALSE)), type = "1", lty = 1, main = "Original data (mean centered)")
matplot(t(scale(Xpv, scale = FALSE)), type = "1", lty = 1, main = "PV-set (mean centered)")</pre>
```

par(oldpar)

pcvpls

#### Procrustes cross-validation for PLS models

#### Description

Procrustes cross-validation for PLS models

#### Usage

```
pcvpls(
   X,
   Y,
   ncomp = min(nrow(X) - 1, ncol(X), 30),
   center = TRUE,
   scale = FALSE,
   cv = list("ven", 4),
   cv.scope = "global"
)
```

#### Arguments

Х	matrix with predictors from the training set.
Υ	vector or matrix with response values from the training set.
ncomp	number of components to use (more than the expected optimal number).
center	logical, to center or not the data sets
scale	logical, to scale or not the data sets
cv	which split method to use for cross-validation (see description for details).
cv.scope	scope for center/scale operations inside CV loop: 'global' — using globally computed mean and std or 'local' — recompute new for each local calibration set.

#### Details

The method computes pseudo-validation matrix Xpv, based on PLS decomposition of calibration set X, y and cross-validation.

Parameter 'cv' defines how to split the rows of the training set. The split is similar to cross-validation splits, as PCV is based on cross-validation. This parameter can have the following values:

1. A number (e.g. 'cv = 4'). In this case this number specifies number of segments for random splits, except 'cv = 1' which is a special case for leave-one-out (full cross-validation).

2. A list with 2 values: 'list("name", nseg)'. In this case '"name"' defines the way to make the split, you can select one of the following: '"loo"' for leave-one-out, '"rand"' for random splits or '"ven"' for Venetian blinds (systematic) splits. The second parameter, 'nseg', is a number of segments for

#### pcvpls

splitting the rows into. For example, 'cv = list("ven", 4)', shown in the code examples above, tells PCV to use Venetian blinds splits with 4 segments.

3. A vector with integer numbers, e.g. 'cv = c(1, 2, 3, 1, 2, 3, 1, 2, 3)'. In this case number of values in this vector must be the same as number of rows in the training set. The values specify which segment a particular row will belong to. In case of the example shown here, it is assumed that you have 9 rows in the calibration set, which will be split into 3 segments. The first segment will consist of measurements from rows 1, 4 and 7.

Parameter 'cv.scope' influences how the Procrustean rule is met. In case of "global" scope, the rule will be met strictly - error of predictions for PV-set and the global model will be identical to the error from conventional cross-validation. In case of "local" scope, every local model will have its own center and hence the rule will be almost met (the errors will be close but not identical).

#### Value

Pseudo-validation matrix (same size as X) with an additional attribute, 'D' which contains the scaling coefficients (ck/c)

#### References

1. S. Kucheryavskiy, O. Rodionova, A. Pomerantsev. Procrustes cross-validation of multivariate regression models. Analytica Chimica Acta, 1255 (2022) [https://doi.org/10.1016/j.aca.2023.341096]

#### Examples

```
# load NIR spectra of Corn samples
data(corn)
X <- corn$spectra
Y <- corn$moisture
# generate Xpv set based on PCA decomposition with A = 20 and venetian blinds split with 4 segments
Xpv <- pcvpls(X, Y, ncomp = 20, center = TRUE, scale = FALSE, cv = list("ven", 4))
# show the original spectra and the PV-set (as is and mean centered)
oldpar <- par(mfrow = c(2, 2))
matplot(t(X), type = "1", lty = 1, main = "Original data")
matplot(t(Xpv), type = "1", lty = 1, main = "PV-set")
matplot(t(scale(X, scale = FALSE)), type = "1", lty = 1, main = "Original data (mean centered)")
matplot(t(scale(Xpv, scale = FALSE)), type = "1", lty = 1, main = "PV-set (mean centered)")
par(oldpar)
```

```
# show the heatmap with the scaling coefficients
plotD(Xpv)
```

PV-

pcvreg

## Description

This is a generic method, use 'pcvpls()' or 'pcvpcr()' instead.

## Usage

```
pcvreg(
  X,
  Y,
  ncomp = min(nrow(X) - 1, ncol(X), 30),
  cv = list("ven", 4),
  center = TRUE,
  scale = FALSE,
  funlist = list(),
  cv.scope = "global"
)
```

## Arguments

Х	matrix with predictors from the training set.
Υ	vector with response values from the training set.
ncomp	number of components to use (more than the expected optimal number).
cv	which split method to use for cross-validation (see description of method 'pcv-pls()' for details).
center	logical, to center or not the data sets
scale	logical, to scale or not the data sets
funlist	list with functions for particular implementation
cv.scope	scope for center/scale operations inside CV loop: 'global' — using globally computed mean and std or 'local' — recompute new for each local calibration set.

plotD	Plots heatmap for scaling coefficients obtained when generating
	set for PCR or PLS

## Description

Plots heatmap for scaling coefficients obtained when generating PV-set for PCR or PLS

## rotationMatrixToX1

## Usage

```
plotD(
   Xpv,
   colmap = colorRampPalette(c("blue", "white", "red"))(256),
   lim = c(-2, 4),
   xlab = "Components",
   ylab = "Segments",
   ...
)
```

#### Arguments

Хрv	PV-set generated by 'pcvpcr()' or 'pcvpls()'.
colmap	colormap - any with 256 colors.
lim	limits for color map (smallest/largest expected value), centered around 1.
xlab	label for x-axis
ylab	label for y-axis
	any other parameters for method 'image'

## Value

No return value, just creates a plot.

rotationMatrixToX1 Creates a rotation matrix to map a vector x to [100...0]

## Description

Creates a rotation matrix to map a vector x to [1 0 0 ... 0]

## Usage

```
rotationMatrixToX1(x)
```

## Arguments

x Vector (sequence with J coordinates)

## Value

Rotation matrix (JxJ)

simpls

## Description

SIMPLS algorithm for calibration of PLS model

## Usage

simpls(X, Y, ncomp)

## Arguments

Х	a matrix with x values (predictors)
Υ	a matrix with y values (responses)
ncomp	number of components to calculate

## Value

a list with computed weights, x- and y-loadings for PLS regression model.

#### References

[1]. S. de Jong. SIMPLS: An Alternative approach to partial least squares regression. Chemometrics and Intelligent Laboratory Systems, 18, 1993 (251-263).

# Index

simpls, 14

15